

## ULT SERIES DIMENSION & WEIGHT

Model	Packing dimensions	Nett/Gross Weight
-------	--------------------	-------------------

### Chest Freezers

DW-86W100J	845 x 855 x 1250	138/160kg
DW-86W420J	2195 x 895 x 1130	310/357kg

### Upright Freezers

DW-86L338J	875 x 970 x 2010	238/278kg
DW-86L338J	875 x 970 x 2010	238/278kg
DW-86L490J	925 x 985 x 2150	295/335kg
DW-86L628	1080 x 965 x 2150	301/323kg
DW-86L578J	950 x 1055 x 2150	300/330kg
DW-86L728J	1090 x 1050 x 2150	345/385kg
DW-86L729		
DW-86L828J	1190 x 1045 x 2150	380/410kg
DW-86L959		

### Ultra Low Energy

DW-86L579BP(T)	950 x 1055 x 2150	325/355kg
DW-86L729BP(T)	1100 x 1105 x 2150	350/385kg
DW-86L829BP(T)	1190 x 1045 x 2150	380/415kg
DW-86L959BP(T)	1365 x 1105 x 2150	450/485kg

### Dual Compressor

DW-86L578S	950 x 1055 x 2150	325/355kg
DW-86L728S	1100 x 1105 x 2150	350/385kg