

Model	Racks: sliding, side load / Chest rack	Rack Model part number	Number of Cryoboxes per rack	Number of racks per shelf	Number of racks per freezer full load
DW-86W100J	Chest Freezer Rack	DCJ-08-A	8 x 2 inch cryobox	N/A	9 racks per freezer
DW-86W420J	Chest Freezer Rack	DCJ-10-A	10 x 2 inch cryobox	N/A	27 racks per freezer
DW-86L338J	Sliding draw	SSDR16-2	16 x 2 inch cryobox	4 racks per shelf	12 x 16 place racks
	Side load	SSSL16-2	16 x 2 inch cryobox	4 racks per shelf	12 x 16 place racks
DW-86L490J	Sliding draw	SSDR16-2	16 x 2 inch cryobox	4 racks per shelf	8 x 16 place + 8 x 20 place racks
	Side load	SSSL16-2	16 x 2 inch cryobox	4 racks per shelf	8 x 16 place + 8 x 20 place racks
	Sliding draw	SSDR20-2	20 x 2 inch cryobox	4 racks per shelf	8 x 16 place + 8 x 20 place racks
	Side load	SSSL20-2	20 x 2 inch cryobox	4 racks per shelf	8 x 16 place + 8 x 20 place racks
DW-86L578J	Sliding draw	SSDR25-2	25 x 2 inch cryobox	4 racks per shelf	16 x 25 place racks
DW-86L579BP	Side load	SSSL25-2	25 x 2 inch cryobox	4 racks per shelf	16 x 25 place racks
DW-86L578S					
DW-86L728J	Sliding draw	SSDR25-2	25 x 2 inch cryobox	5 racks per shelf	20 x 25 place racks
DW-86L729BP	Side load	SSSL25-2	25 x 2 inch cryobox	5 racks per shelf	20 x 25 place racks
DW-86L728S					
DW-86L828J	Sliding draw	SSDR25-2	25 x 2 inch cryobox	6 racks per shelf	24 x 25 place racks
DW-86L829BP	Side load	SSSL25-2	25 x 2 inch cryobox	6 racks per shelf	24 x 25 place racks
DW-86L959BP	Sliding draw	SSDR25-2	25 x 2 inch cryobox	7 racks per shelf	28 x 25 place racks
	Side load	SSSL25-2	25 x 2 inch cryobox	7 racks per shelf	28 x 25 place racks