

Pharmacy Refrigerators

— HYC-68/68A



Ergonomic Design



Audible and Visual Alarms

Audible buzzer and visible flashing light



Multiple Alarm Functions

High/low temperature, power failure, door ajar, sensor error, low battery, and remote alarm



Air Cooling, Auto Defrost

Superior cabinet temperature uniformity, auto defrost



Accurate Temperature Control



Interior LED



Lockable Design

Qingdao Haier Biomedical Co.,Ltd.

No.280 Feng Yuan Road, High-tech Zone,
Qingdao, 266109, P.R. China
Tel: +86-0532-88935593
Website: www.haiermedical.com



Haier Biomedical
International



Haier Biomedical
International



@haiermedicalint



Haier Biomedical
International



Haier Biomedical
International

Specifications

Technical Data

Cabinet Type	Upright
Climate Class	N
Cooling Type	Forced air cooling
Defrost Mode	Auto
Refrigerant	HC
Noise (dB (A))	41

Performance

Temp Range (°C)	2~8
-----------------	-----

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply (V/Hz)	220~240/50	115/60
Power (W)	140	150
Electrical Current (A)	0.85	1.7

Dimensions

Capacity (L/Cu.Ft)	68/2.4
Net/Gross Weight (approx)	38/40 (kg) 83.8/88.2 (lbs)
Interior Dimensions (W*D*H)	415*385*505 (mm) 16.3*15.2*19.9 (in)
Exterior Dimensions (W*D*H)	495*580*660 (mm) 19.5*22.8*26.0 (in)
Packing Dimensions (W*D*H)	540*640*710 (mm) 21.3*25.2*28.0 (in)
Container Load (20'/40'/40'H)	102/204/204



HYC-68



HYC-68A

Functions

High/Low Temp	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
Door Ajar	Y
Remote Alarm	Y

Accessories

Caster	Y
Foot	Y
Porthole	Y
Shelves /Drawers	2/1

Other

Certification	CE	UL
---------------	----	----